

Case Competition

Eligibility: General Practitioner • General Cardiologist
Interventional Cardiologist

Guideline

- 1 The case report should be presented in structured abstract as below:

Background

Case Illustration
/ Presentation*

Conclusion

Keyword
(no more than 5)

*Discussion is not required; however, if needed can be merged into the 'Case illustration' section with the same word limit of 350.

- 2 The case report should not be exceeded 350 words (*excluding the title, authors, and affiliations*).
- 3 The authors' names must adhere to the following format: first & middle name written as initials, last name in full (e.g., *R.P. Chandra*) with affiliations.
- 4 Maximum of 3 figures (1 is mandatory) OR table is allowed to be included in the articles (*in the form of PNG, JPEG, or JPG. max 15MB*).
- 5 Maximum of 2 videos are allowed to be included in the articles with the maximum size of 100MB per video.

Submit before 14 Sept 2024

For more info please contact: +62 813-1522-7958 (Santi) or follow our IG @isicam_inalive

Research Competition **META-ANALYSIS** **/SYSTEMATIC REVIEW**

Eligibility: General Practitioner • General Cardiologist
Interventional Cardiologist

Guideline

- 1 The case report should be presented in structured abstract as below:

Background	Method	Result	Conclusion	Keyword (no more than 5)
------------	--------	--------	------------	-----------------------------
- 2 The case report should not be exceeded 350 words (*excluding the title, authors, and affiliations*).
- 3 The authors' names must adhere to the following format: first & middle name written as initials, last name in full (e.g., *R.P. Chandra*) with affiliations.
- 4 A mandatory of 1 figure (max 15MB) and 1 .docx file (max 15MB) is needed to be included in the articles (*in the form of PNG, JPEG, or JPG and docx*).
- 5 Maximum of 2 videos are allowed to be included in the articles with the maximum size of 100MB per video.

Submit before 14 Sept 2024

For more info please contact: +62 813-1522-7958 (Santi) or follow our IG @isicam_inalive